
# Symphony®

EN Instructions for use | ES Instrucciones de uso | FR Mode d'emploi

Read all instructions before using this product.

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# Table of Contents

<b>1. Important Safeguards</b> .....	<b>4-5</b>
<b>2. Product Description</b> .....	<b>6-9</b>
<b>3. Getting Started</b> .....	<b>10</b>
<b>4. Cleaning Overview</b> .....	<b>11</b>
<b>5. Wash - Before first use &amp; after each use</b> .....	<b>12-13</b>
<b>6. Sanitize - Before first use &amp; once per day</b> .....	<b>14-15</b>
<b>7. Tubing and Breastpump Care</b> .....	<b>16-17</b>
7A. Tubing Care.....	16-17
7B. Cleaning and Disinfecting the Breastpump.....	17
<b>8. Powering Your Pump</b> .....	<b>18-19</b>
8A. Card and Cord Protector.....	18-19
8B. Powering On.....	18-19
8C. Changing the Program Card.....	19
8D. Traveling Outside the U.S.....	19
<b>9. Assembly of the Pump Kit</b> .....	<b>20-23</b>
<b>10. Breastshield Sizing</b> .....	<b>24</b>
<b>11. Pumping</b> .....	<b>25-27</b>
11A. Before pumping.....	25
11B. Preparing for pumping.....	25-27
11C. Finish pumping.....	26-27
<b>12. Premie+™ Initiation Card</b> .....	<b>28-29</b>
<b>13. Storing, Preparing &amp; Feeding Breastmilk</b> .....	<b>30</b>
13A. Storing Breastmilk.....	30
13B. Preparing Breastmilk.....	30
13C. Feeding Breastmilk.....	30
<b>14. Troubleshooting</b> .....	<b>31-32</b>
<b>15. Error Codes</b> .....	<b>33</b>
<b>16. Supplemental Information</b> .....	<b>34</b>
<b>17. Warranty</b> .....	<b>35</b>
<b>18. Meaning of Symbols</b> .....	<b>36-37</b>
<b>19. EMC/Technical Description</b> .....	<b>38-43</b>
<b>20. Technical Specifications</b> .....	<b>44</b>
<b>21. Disposal</b> .....	<b>45</b>

## IMPORTANT

- Plastic bottles and component parts become brittle when frozen and may break when dropped.
- Bottles and component parts may become damaged if mishandled, e.g. dropped, over-tightened, or knocked over.
- Take appropriate care in handling bottles and components.
- Do not use the breastmilk if bottles or components become damaged.

## Intended Use

The Symphony® Breastpump is a powered breastpump to be used by lactating women to express and collect milk from their breast.

## Product Description

Symphony is a multi-user, hospital grade electric breastpump that includes 2-Phase Expression® technology and is capable of single and double pumping.

The warning symbol identifies all instructions that are important to safety. Failure to observe these instructions can lead to injury or damage to the breastpump. When used in conjunction with the following words, the warning symbols stand for:

**WARNING** Can lead to serious injury or death.

**CAUTION** Can lead to minor injury.

**NOTE** Can lead to material damage.

**TIP** Useful or important information that is not related to safety.

**PLEASE SAVE THESE INSTRUCTIONS.**

# 1. Important Safeguards

When using electrical devices, basic safety precautions should always be followed, including the following:

## **READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.**

### **DANGER: To reduce the risk of electrocution:**

- Always unplug electrical product immediately after use.
- Do not use while bathing or showering.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not submerge pump unit into water or other liquid.
- Do not reach for a product if it has fallen into water or other liquid. Immediately unplug from electrical outlet and contact manufacturer.
- Do not spray or pour liquid directly onto the pump unit.

### **⚠ WARNING: To avoid fire, electrocution, or serious burns:**

- Do not leave product unattended when plugged into an electrical outlet.
- Close supervision is necessary when this product is used near children.
- Use the product only for its intended use as described in this manual.
- Do not use attachments not recommended by the manufacturer, as this may cause hazards.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Keep power cord and all components away from heated surfaces.
- Never use while sleeping or when drowsy.
- Never drop or insert any object into any opening or tubing.
- Do not use outdoors.
- Keep appliance dry and do not place in direct sunlight.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Separation from the mains is only assured through the disconnection of the power cord from the wall socket outlet.
- Position the device to ensure easy disconnection of the power cord from the supply mains.
- Use Medela original accessories only.
- There are no user-serviceable parts inside the Symphony. DO NOT ATTEMPT TO REPAIR THE Symphony Breastpump YOURSELF. No modifications to the device are permitted.
- Always inspect power cord wires prior to use for damage or exposed wire. If damage is found, immediately discontinue use of power cord and call Medela Customer Service at 1-800-435-8316.

### **WARNING: To avoid health risk and reduce the risk of injury:**


- Breastpump kits are personal care items intended for use by a single user only and should not be shared between users.
- Inspect the breastpump kit components prior to each use.
  - Do not use the components if they are damaged.
  - Do not use the components if dirt, mold or other contamination is visible.
- Inspect all appropriate pump components before each use.
- Clean and sanitize all parts that come in contact with your breast and breastmilk prior to first use.
- Do not modify any breastpump kit components.
- Do not drive while pumping.
- Never use while pregnant, as pumping can induce labor.
- Wash all parts that come in contact with your breast and breastmilk after every use.
- Do NOT continue pumping for more than 2 consecutive pumping sessions if no results are achieved.
- Do not thaw frozen breastmilk in a microwave or in a pan of boiling water.
- If tubing becomes moldy, discontinue use and replace tubing.
- If breastmilk backup occurs you must clean tubing before your next pumping session.
- If problems or pain occurs, consult your lactation consultant or physician.

### **CAUTION: Can lead to minor injury:**

- Pump only with the breastpump in an upright position.
- Do NOT shorten tubing between breastshield and pump.
- Make sure the voltage of the device is compatible with the power source.
- Plug the power cord into the breastpump first and then into the wall socket.
- Do NOT use abrasive cleaners/detergents when cleaning breastpump or breastpump parts.
- Never put breastpump in water or in a sterilizer, as you can cause permanent damage to the breastpump.
- Do not attempt to remove the breastshield from your breast while pumping. Turn the breastpump off and break the seal between your breast and breastshield with your finger, then remove breastshield from your breast.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between the breast and the breastshield with your finger and remove the breastshield from your breast.
- Do not store wet or damp parts as mold may develop.
- Do not use the Medela Symphony Breastpump kit with a non-compatible breastpump.
- Do not use the breastpump kit for any other purpose than expressing breastmilk.
- Portable and mobile radio frequency communications equipment can affect the Symphony Breastpump.
- The Symphony Breastpump needs special precautions regarding EMC (Electromagnetic compatibility) and needs to be installed and put into service according to the EMC information provided in these instructions for use.

## 2. Product Description

Authentic Medela Spare Parts are made with the quality materials and workmanship you expect and trust, only from Medela.

Non-Medela spare parts can vary significantly in terms of product design, materials, and workmanship; all of which may affect the performance of your Medela breastpump.

Quality & performance  
you can trust

## Symphony Breastpump

0240110  
0240108  
0240208

On/Off button

Vacuum adjustment

LCD-display

Let-down button

medela


Container Stand  
8100552

Card & Cord  
Protector  
8100589  
Replacement Screw  
4130003

Power cord  
9280043

Socket for 12Volt  
Vehicle Lighter  
Adaptor (sold separately)  
67173

Program card  
Standard 2.0

67208 - English  
(Included with 0240108 and 0240208)

67209 - French  
67210 - Spanish

Program card Premie<sup>+</sup>

67227 - English  
(Included with 0240110)

27225 - French  
67226 - Spanish

French and Spanish Program Cards  
(sold separately)

Case  
6007080  
(sold separately)

## 2. Product Description

### Optional Accessories

Symphony® Hard Case	6007080
Symphony Trolley	6007092
Symphony Vehicle Lighter Adaptor	67173
Symphony Cooler Carrier - Black	37002
Breastmilk Transport Bag™ - NICU	67235-06
Cavicide®	3007009
Breastmilk Pumping Log - English	1547553
Breastmilk Pumping Log - Spanish	1547556

### Breastpump Kits


Symphony Double Pumping Kit	67099-06
Symphony & Lactina® Double Pumping Kit - Sterile	67316S
Symphony Double Pumping Kit - Sterile	67399S
Symphony Breastmilk Initiation Kit - Sterile	67350S
Symphony Harmony™ Double Duet Breastpump Kit - Sterile	67115S
Symphony Harmony Double Duet Breastmilk Initiation Kit - Sterile	67355S
Symphony Lactina Breastmilk Initiation Kit - Sterile	67340S
Symphony Single Pumping Kit - Sterile	67026S
Symphony Harmony Duet Breastpump Kit - Sterile	67114S
Symphony Lactina Double Pumping Kit - Retail	67116-06
Symphony & Lactina Rental Kit	67206-03
Symphony to Lactina Conversion Kit	67090
Lactina to Symphony Conversion Kit	67091
Harmony to Symphony Conversion Kit	67023

## Replacement Parts

Symphony® Diaphragm Kit	8007199
Symphony Container Stand	8100552
Symphony Membrane Cap	8100769
Symphony Protective Membrane	8200094
Symphony Power Cord	9280043
Symphony Preemie+™ Card English	67227
Symphony Preemie+ Card Spanish	67226
Symphony Preemie+ Card French	27225
Symphony Program Card/Power Cord Protector	8100589
Symphony Program Card/Power Cord Protector Replacement Screw	4130003
Symphony 2.0 Program Card - English	67208
Symphony 2.0 Program Card - French	67209
Symphony 2.0 Program Card - Spanish	67210
PersonalFit™ Connector	87071
PersonalFit 21mm Breastshield	87072
PersonalFit 24mm Breastshield	87073
PersonalFit 27mm Breastshield	87274
PersonalFit 30mm Breastshield	87075
PersonalFit 36mm Breastshield	87084
XL Glass Breastshield - 40mm	6107041
Replacement Membranes	87088
Replacement Valves & Membranes	87089

## 3. Getting Started

### Note

It's important for you to do the following before using Symphony® for the first time:

1. For models that include a rechargeable battery, fully charge the battery 12 hours uninterrupted. For further battery care information see **Section 8**.
2. Separate all parts
3. Clean – see instructions:
  - Wash – **Section 5**, pages 12-13
  - Sanitize – **Section 6**, pages 14-15

Breastshield assemblies come assembled. Separate all parts that will come into contact with the breast and breastmilk before cleaning.

#### Parts to clean:

- Breastshields
- Breastmilk bottles
- Lids
- Connectors
- Valves
- Membranes

## 4. Cleaning Overview

	Breastpump kit & bottles 	Tubing 	Symphony® Breastpump and Case 
<b>Before 1st use</b>	<ul style="list-style-type: none"> <li>✓ wash</li> <li>✓ sanitize</li> </ul>		
<b>After each use</b>	<ul style="list-style-type: none"> <li>✓ wash</li> </ul>		
<b>Once per day</b>	<ul style="list-style-type: none"> <li>✓ sanitize</li> </ul>		
<b>As needed</b>		<ul style="list-style-type: none"> <li>✓ wash or sanitize <b>only if</b> milk or condensation in tubing</li> </ul>	<ul style="list-style-type: none"> <li>✓ wipe with clean, damp cloth</li> </ul>

Detailed cleaning instructions on pages 12-17:

- Wash – page 12-13
- Sanitize – page 14-15
- Tubing and breastpump care – pages 16-17

### Note

Sterile kits are ready to use and do not require washing prior to first use.

### Warning

Breastpump kits are a personal care item intended for use by a single user only and should not be shared between users.

## 5. Wash – Before first use and after each use\*

### Supplies needed:

- Mild dish soap
- Clean dish towel or soft brush
- Clean sink or bowl
- Drinking-quality water

### Parts to wash:

- Breastshields
- Breastmilk bottles
- Lids
- Connectors
- Valves
- Membranes

### ! Note

\*Sterile kits are ready to use and do not require washing prior to first use.

### ! Warning

Breastpump kits are a personal care item intended for use by a single user only and should not be shared between users.

1

- Separate parts for washing

### ! Caution

Separate and wash all parts that come in contact with breastmilk immediately after use to avoid dry up of breastmilk residue and prevent growth of bacteria.

### ! Note

- Wash breastpump kit parts after each use.
- Take care not to damage parts of the pump kit during cleaning.
- **It is not necessary to wash tubing, membrane cap or protective membrane prior to first use.**
- **Only wash tubing if condensation or milk is present.**

## Wash in dishwasher

4

- Wash all separated parts on top rack of dishwasher.
- Allow all breastpump parts to air dry in a clean area.
- Store dry parts in a clean, cool place when not in use.

### ! Note

If the individual components of the pump kit are cleaned in a dishwasher, parts may be discolored by food pigments. This will not impact part function.

## Wash in sink

2

Rinse all separated parts that came in contact with breast and breastmilk in cool water in order to remove breastmilk residue.

3

OR

- Soak all separated parts in warm, soapy water for 5 minutes.
- Wash each part with a clean dish-cloth or soft brush.
- Rinse all separated parts with clear water.

### Caution

Only use drinking-quality tap or bottled water for cleaning.

*You may also wash your kit parts by following the instructions on the Quick Clean™ Breastmilk Removal Soap bottle. (sold separately)*

5

- Place parts on a clean surface and/or towel.
- Allow all parts to air dry.
- Store dry parts in a clean, cool place when not in use. Do NOT store wet or damp parts.

### Note

Store the pump kit in a clean bag/container until next use.

## 6. Sanitize – Before first use and once per day\*

### Supplies needed:

- Clean pot for boiling water
- Drinking-quality water
- Tongs

### Parts to sanitize:

- Breastshields
- Breastmilk bottles
- Lids
- Connectors
- Valves
- Membranes

1

- Separate all parts that come in contact with breastmilk.

### Note

Sanitize breastpump kit parts once daily.

### Note

\*Sterile kits are ready to use and do not require washing or sanitizing prior to first use.

### Warning

Breastpump kits are a personal care item intended for use by a single user only and should not be shared between users.

4

Allow water to cool and gently remove parts from water with tongs.

- Wash hands thoroughly.

- Fill a pot with enough water to cover all parts.
- Bring water to a boil.
- Place parts in boiling water for 10 minutes.

**!** Note

If you notice a white residue on your parts after boiling, you may have a high mineral content in your water. Remove residue by wiping parts with a clean washcloth and air dry. We recommend that you boil parts with distilled water to prevent substantial mineral build-up over time which may compromise your parts.

- Place parts on a clean surface and/or towel.
- Allow all parts to air dry.
- Store dry parts in a clean, cool place when not in use. Do NOT store wet or damp parts.

*You may also sanitize your kit parts by following the instructions on Medela's Quick Clean™ Micro-Steam™ bag. (sold separately)*

**!** Note

Store the pump kit in a clean bag/container until next use.

## 7. Tubing & Breastpump Care

### Supplies needed:

- Mild dish soap
- Clean dish towel
- Clean sink or bowl
- Drinking quality water
- Damp (not wet cloth)

### Parts needed:

- Tubing
- Symphony® pump
- Membrane cap
- Protective membrane

### 7A. Tubing care

- Inspect tubing after each pumping session.

#### Caution

- Do not store wet or damp parts as mold may develop.
- If tubing becomes moldy, discontinue use and replace tubing. Contact Medela Customer Service at customer.service@medela.com or call 1-800-435-8316.

#### Note

**Tubing should only be washed if there is any sign of condensation or milk.**

3

Use a damp (not wet) cloth to wipe the pump and membrane area.

4

Remove, tubing, protective membrane, and membrane cap. Clean following steps in section 5.

#### Warning

Breastpump kits are a personal care item intended for use by a single user only and should not be shared between users.

7


- Wash tubing in warm soapy water.
- Rinse tubing with clear water.
- Hang to air dry.
- Or shake out any drops of water from the tubing and let the pump run with tubing attached for 1-2 minutes until dry.

1 Turn off breastpump.

2 Unplug breastpump from power source.

5 Remove tubing from breastshield.

6 Rinse tubing in cool water to remove breastmilk.

## 7B. Cleaning and disinfecting the breastpump

1 Wipe off with a clean, damp (not wet) cloth. Water and basic dishwasher agent or non-abrasive detergent may be used.

Wipe off with a clean, damp (not wet) cloth. Alcohol based disinfecting agent may be used. Do not immerse the pump in water; do not run water over the pump.

### Caution

Unplug breastpump before cleaning.

## 8. Powering Your Pump

### Note

- Only use the pump with the appropriate program card (Standard or Preemie+).

### Optional for pump with battery:

- Built-in battery must be charged for 12 hours on receipt of the device.
- The pump can be used during charging.
- If the pump is to remain unused, charge the battery every two months (for 12 hours each time).
- For optimum service life of the battery, charge the battery while breastpump is not in use and avoid fully discharging the battery on a regular basis.

Pumping time (fully-charged battery)	approx. 60 min.
Charging time (for full charge)	12 hours

### 8A. Card and Cord Protector


To remove, unscrew the card and cord protector and slide it out of the guide rail.

### 8B. Powering On

Insert the power cord into the power socket on the back of the pump.


Insert the cord into the card and cord protector in a loop around the pin.

→ Make sure there is enough play to plug the cord in.

Push the card and cord protector over the guide rail on the back until it stops and tighten the screw.

### 8C. Changing the program card

Plug the power cord into the mains socket.

As soon as the card compartment is revealed, the card can be removed and replaced or changed.

#### Warning

- Only use the power cord that comes with the Symphony®.
- Make sure the voltage of the power cord is compatible with the power source.

### 8D. Traveling outside of the U.S.

While travelling internationally, we recommend purchasing a Universal Power Plug Adaptor (not included). Please consult with the country you are visiting to find out what type of adaptor will work best. Symphony power cord is two-pronged and does not require a grounded outlet.

For information about traveling with your breastpump and expressed breastmilk, visit [www.TSA.gov](http://www.TSA.gov).

## 9. Assembly of Pump Kit

### **T** Tip

Carry out all the steps with care and assemble the pump kit correctly for optimum performance.

### Parts Needed:

- Tubing
- Symphony® pump
- Breastshields
- Breastmilk bottles
- Lids
- Connectors
- Valves
- Membranes
- Membrane cap
- Protective membrane

### **!** Note

- Use Authentic Medela Spare Parts only. See page 6 for details.
- Check pump kit components for wear or damage before use and replace if necessary.
- All components must be completely dry before use for proper performance and to prevent damage to the pump.

### **!** Caution


Wash hands thoroughly with soap and water before touching breastpump, kit and breasts and avoid touching the inside of containers or lids.

Push the valve head and membrane onto the connector.  
→ Position valve head laterally.

Place the protective membrane onto a flat surface and carefully push the membrane cap down until it engages.

2  
Carefully push the membrane onto the valve head until it engages.


3  
Push the breastshield onto the connector.

5  
Screw the bottle into the connector.

6  
Insert the tubing into the corresponding opening in the connector.


8  
Insert the tubing adapter into the opening in the membrane cap.

## 9. Assembly of Pump Kit (cont.)

- Open the lid of the breast pump.  
→ Press the oval button on the top of the Symphony and lift the handle.

- Position the membrane cap firmly on the breast pump membrane to form a seal. Center the tubing connector in the recess in the pump housing.  
→ Make sure that the membrane cap engages and lies flush to the surface.

Center the tubing connector in the recess in the pump housing.  
→ Make sure that the membrane cap engages and lies flush to the surface.

Close the lid.

#### Note

- **For step 12:**

Always close the lid when pumping. The lid keeps the membrane cap pushed down to prevent a loss of vacuum that could occur if the membrane cap came off.

# 10. Breastshield Sizing

## Caution

- While some discomfort may be felt when first using a breastpump, using a breastpump should not cause pain.
- If you are experiencing discomfort at the base of the nipple due to rubbing of your breast tissue with the breastshield tunnel, use of a lubricant such as Tender Care Lanolin may be beneficial.
- For assistance with correct breastshield sizing and comfort please visit [www.medelabreastshields.com](http://www.medelabreastshields.com) or see a lactation consultant / breastfeeding specialist.

## Correctly fitting breastshields

Properly fitting breastshields will ensure comfort and efficient milk expression. To determine if you have the correct size, look at the following diagram:

- The nipple should be centered in the tunnel of the shield.
- The nipple should move freely without rubbing when the pump is turned on.
- Minimal or no part of your areola should be pulled into the tunnel of the breastshield.
- You should see a gentle, rhythmic motion in the breast with each cycle of the pump.
- After pumping, your breast should feel much less full with no areas of hardness.

**Correct fit**

## PersonalFit™ Breastshields

**21 mm**  
Item 87072

**24 mm** *Included with  
Item 87073 *Symphony® Kits**

**27 mm**  
Item 87274

**30 mm**  
Item 87075

**36 mm**  
Item 87084

### Choosing the right size PersonalFit™ breastshield:


1. Determine the size you're currently using. If you're not sure, look for the size embossed on your breastshield (see picture). 24 mm breastshields are provided with Medela breastpumps.

2. While pumping, compare your fit to the images below as a sizing guideline.

**Correct fit**

Your breastshield is **too small**;  
try a larger size

Your breastshield is **too large**;  
try a smaller size

**PersonalFit™ Breastshields are available at many local retailers or visit [www.shopmedela.com](http://www.shopmedela.com).**

# 11. Pumping

## ! Note

- If the pump has been running for 30 minutes without any user action (e.g. adjustment of the vacuum), the Symphony will automatically switch off.
- Only fill the bottle to the 150 mL mark.

## T Tip

- Do not tilt or overfill containers when pumping.


## 11A. Before pumping

## ⚠ Caution

Wash hands thoroughly with soap and water before touching breastpump, kit and breasts and avoid touching the inside of containers or lids.

## 11B. Preparing for pumping

1 Push membrane caps into place on the pump. A "click" will be heard once properly seated. (**Note:** if single pumping, only seat one membrane cap.)

2 Place the breastshield on your breast so that your nipple is properly **centered** in the tunnel. Hold the breastshield onto your breast with your thumb and index finger. Support your breast with the palm of your hand.

## ⚠ Caution

Pressing shields too firmly onto breasts may affect milk flow.

## 11. Pumping (cont.)

Switch on the breastpump with .  
 → Stimulation phase begins.  
 The stimulation phase will automatically change to the expression phase after two minutes.

Place the second breastshield on the breast so that the nipple is properly **centered** in the tunnel.

### 11C. Finish Pumping

Switch off the breastpump with .  
 → Do not fill the bottle above maximum volume indicator line. Always disconnect the Symphony® from the power source after expression.

Use the container stand to prevent the bottle from tipping over.

### Glossary of Terms

<b>2-Phase Expression®</b> 	Through research, Medela found that babies nurse in 2 Phases – Stimulation and Expression. This research is the basis for the technology in all Medela breastpumps.
<b>Stimulation Phase</b>	Fast sucking/pumping rhythm to stimulate milk flow.
<b>Expression Phase</b>	Slower sucking/pumping rhythm for gentle and efficient milk removal after milk has started flowing.
<b>Maximum Comfort Vacuum™</b>	Highest vacuum setting where pumping still feels comfortable. Different for every mother.

If your milk starts to flow before the end of the two minutes, press  to change to the expression phase.  
 → Make sure that the milk flows into the bottle properly.

Set vacuum comfort: Rotate the vacuum to the right and increase until there is slight discomfort. Rotate to the left to reduce the vacuum slightly.

Close the bottle with a lid (depending on bottle type).  
 → Follow the instructions in

**Section 13** “Storing, preparing and feeding breastmilk”.

**If you are experiencing low or no suction, refer to Troubleshooting Section 14, pages 31-32.**

### Caution

- Contact your healthcare professional or lactation professional if you can express only minimal or no milk or if expression is painful. See **Section 10 and 16** for more information.
- Do not try and express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.

### Note

- Press the “Let down” button  to revert from the expression phase to the stimulation phase if necessary. If milk has not begun to flow after two successive stimulation phases, take a break and try again in 15–30 minutes. Try breast massage (ask your lactation professional or health care practitioners) and then pump using the expression phase.
- The program stores the last vacuum setting in the stimulation phase.


## 12. Premie+™ Initiation Card

### Note

- Assists in initiating and maintaining lactation in breastpump-dependent mothers of premature infants.
- The Premie+ Initiation Card contains both pumping programs.
- **Premie+ 1.0:**  
Program for the immediate post-natal period until initiation of milk production. Then use the Standard program. The Premie+ program runs uninterrupted for 15 minutes and comprises a sequence of various stimulation, expression and pause phases. “Done” appears in the display when the program is complete.
- **Standard 2.0:**  
Use after initiation of milk production. Initiation is reached once three consecutive pumping sessions totalling 20 mL of breastmilk (from both breasts combined) is reached. Includes “2-Phase Expression”, starting with the stimulation phase and followed by the expression phase.

1  
Unscrew the card and cord protector (optional) and slide out of the guide rail.

4  
Premie+: switch on the breastpump with . Then press  within 10 seconds.

2

Insert the Preemie+ Initiation card.

3

Label the pump with the Preemie+ sticker. Replace the card and cord protector (optional).

5

Set vacuum comfort: Rotate the vacuum to the right and increase until there is slight discomfort. Rotate to the left to reduce the vacuum slightly.

6

Use the breastmilk pumping log to check and track expression sessions. (optional and sold separately)

**For more details regarding pumping see Section 11.**

## 13. Storing, Preparing & Feeding Breastmilk

### 13A. Storing breastmilk

Freshly Expressed Breastmilk Storage Guidelines (For Healthy Term Babies)			
Room Temperature	Cooler with Ice Pack	Refrigerator	Freezer
4–6 hours at 66–78 °F (19–26 °C)	24 hours at 59 °F (15 °C)	3–8 days at 39 °F or lower (4 °C)	6–12 months 0–4 °F (-18–-20 °C)

References: [www.BreastmilkGuidelines.com](http://www.BreastmilkGuidelines.com)

### 13B. Preparing breastmilk

#### Caution

- Do not thaw frozen breastmilk in a microwave or in a pan of boiling water.
- Do not microwave breastmilk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breastmilk.)
- Thaw breastmilk overnight in the refrigerator. Thawed breastmilk is safe in the refrigerator for 24 hours. Do not refreeze thawed breastmilk.
- Quickly thaw breastmilk by holding the bottle under warm running water.
- Place the sealed bottle in a bowl of warm water for 20 minutes to bring it to body temperature.
- If you are adding expressed breastmilk to a container of already frozen breastmilk, make sure to add a lesser amount than the already frozen amount.

### 13C. Feeding breastmilk

It is recommended that breastfeeding is well established prior to bottle feeding your baby.

- Always inspect the bottle, nipple and other components immediately before and after each use. If nipple appears cracked or torn, discontinue use immediately.
- To prevent possible choking hazard, test strength of nipple by pulling on bulb portion of the nipple.
- Do not attempt to enlarge the nipple hole.
- Infants must not be bottle fed without adult supervision.
- Nipple should not be used as a pacifier.

## 14. Troubleshooting

<p><b>If the breastpump is not running</b></p>	<p>Check whether electrical outlet power/battery is available.</p> <p>Check that the device is switched on.</p> <p>Check whether program card is inserted correctly in the card compartment on the back of the breastpump.</p>
<p><b>If there is low or no suction</b></p>	<p>Assemble the pump set in accordance with the instructions in <b>Section 9</b>.</p> <p>Check whether all connections of the pump set and breastpump are secure.</p> <p>Pump in accordance with <b>Section 11</b>. Make sure that the breastshield forms a complete seal around the breast.</p> <p>The white valve membrane must lie flat against the yellow valve head.</p> <p>Check that the valve head and membrane are clean and not damaged.</p> <p>See <b>Section 9</b> "Assembly of Pump Kit".</p> <ul style="list-style-type: none"> <li>• The ends of the tubing must be firmly connected to the connector and the opening in the membrane cap.</li> <li>• The protective membrane and membrane cap should be correctly assembled and the membrane cap must be firmly fitted to the breastpump membrane to form a seal.</li> <li>• The protective membrane must be intact.</li> <li>• Do not kink or pinch the tubing during expression.</li> <li>• All components must be completely dry.</li> </ul> <p><b>Valve and Membrane Check</b></p> <ul style="list-style-type: none"> <li>• Remove and separate valve from the membrane.</li> <li>• Inspect the valve for chips or cracks. Inspect the membrane for holes or tears. If flawed or damaged, replace the valve and/or membrane prior to pumping.</li> <li>• Wash valve and membrane per Cleaning Instructions in <b>Section 5</b>.</li> <li>• Reassemble. Make sure the membrane lays flat on the valve, and if it does not, replace the valve and/or membrane prior to pumping.</li> </ul> <p><b>Breastshield and Breastshield Connector Check</b></p> <ul style="list-style-type: none"> <li>• Inspect breastshield and breastshield connector to make sure there are no cracks and they securely attach to each other.</li> <li>• Make sure the breastshield forms a complete seal around the breast.</li> </ul> <p><b>Other Checks</b></p> <p>Inspect the following to check for damage and secure connections.</p> <ul style="list-style-type: none"> <li>• Secure tubing adapter to the breastshield connector.</li> <li>• Test vacuum. If suction is not improved, contact Medela Customer Service at 1-800-435-8316.</li> </ul> <p><b>Tubing Check</b></p> <ul style="list-style-type: none"> <li>• If damaged replace prior to pumping.</li> <li>• Inspect tubing for damage and secure connections.</li> <li>• Check that the ends of the tubing are firmly connected to the connector and the opening in the membrane cap.</li> <li>• If single pumping, make sure only one membrane cap is seated.</li> </ul> <p><b>Take care not to kink or pinch tubing.</b></p>

# 14. Troubleshooting

<b>If there is an overflow</b>	<ol style="list-style-type: none"><li>1. Place the pump set into the bottle holder.</li><li>2. Turn off the pump and unplug the power cord from the electrical outlet.</li><li>3. Use a damp (not wet) cloth to wipe the pump and membrane area.</li><li>4. If there is milk in the tubing: Remove tubing, protective membrane and membrane cap and clean in accordance with <b>Section 7</b>. Shake out any drops of water from the tubing and hang up the tubing to dry or let the pump run with tubing attached for 1–2 minutes or until dry.  If there is condensation in the tubing: Let the pump run with tubing attached for 1–2 minutes or until dry.</li><li>5. When components are completely dry, reassemble breastpump and accessories set.</li></ol>
<b>Condensation or Milk in Tubing</b>	<ul style="list-style-type: none"><li>• Turn off the pump and unplug from power source.</li><li>• Use a damp (not wet) cloth to wipe the pump.</li><li>• Refer to cleaning instructions "Tubing Cleaning" <b>Section 7</b>.</li><li>• If necessary, clean the pump kit.</li><li>• When components are completely dry, reassemble breastpump and kit.</li></ul>
<b>No "Let-Down" or Breastmilk Expression</b>	<ul style="list-style-type: none"><li>• Ensure your breastpump kit is assembled correctly and suction is available.</li><li>• Relax and take a 10-15 minute break, if let down is not achieved after 2 consecutive pumping sessions.</li><li>• Consult with a healthcare or lactation professional if expression does not occur.</li></ul>

If you have not resolved the problem with your breastpump or you have further questions, please contact Medela Customer Service at 1-800-435-8316 or fill out the contact us form at [www.medelabreastfeedingus.com](http://www.medelabreastfeedingus.com).

# 15. Error Codes

There are two ways in which a Symphony® Breastpump will communicate that the unit has an issue. If the Symphony is experiencing an issue, an error code will be displayed on the LCD screen (see fig. 1). Previous problems are recorded in the Menu function under “Err:”. To retrieve a history of errors, hold down the “On” button while turning the vacuum control knob, then continue to move the knob until you reach “Err” (see fig. 2). If you continue to experience issues with your Symphony, please contact Customer Service at 1-800-435-8316.

LCD Display Messages		
The following describes LCD messages that may be displayed indicating an issue with the unit, the potential causes, and action steps to correct the issue.		
	Potential Issue	Actions
Invalid card	Program card isn't a Symphony card	Check if correct program card is used
	Program card is incorrectly inserted	Check that card is positioned correctly
	Program card is damaged and not working	Replace program card
No card	Program card isn't inserted	Insert program card
	Program card isn't completely inserted	Push card in all the way
Battery low	Battery needs charging soon, will beep every 20 seconds (15 minutes of pumping time remaining from when the indication started)	Be sure power cord is plugged into back of pump
Battery flat	Battery is completely discharged	Charge the batteries, but you may continue to pump
Motor overload	Motor is drawing too much current	Send Symphony pump to Medela® service center if this happens repeatedly
Configuration	Software configuration	Replace program card if experiencing error (See Technical Bulletin No. 109A)
Error	Potential fault in electronics of gear control	Send Symphony pump to Medela service center if this happens repeatedly
Power fault	Potential fault in electronics of power supply	Send Symphony pump to Medela service center if this happens repeatedly

## 16. Supplemental Information

### Supplemental Information

Following are some common breastfeeding related conditions. If you experience any of these symptoms, contact a healthcare professional or lactation professional.

	Symptom	Potential Causes
<b>Engorgement</b>	Breasts are hard and uncomfortable, possibly with reddened areas.	Milk not draining properly from breast. Milk can drain back into the tissue causing areas to swell and become tender.
<b>Blocked/Plugged Ducts</b>	Area of breast will look reddened and may be tender to touch.	Milk not being drained from a specific duct. The area becomes "clogged" up and milk is then prevented from flowing.
<b>Mastitis</b>	It often follows engorgement but can just occur suddenly. The initial symptoms are similar to flu – fatigue, headache and muscle aches, fever and localized breast tenderness. It can be in one or even both breasts and needs immediate treatment.	A bacterial infection in the breast tissue which is often accompanied by cracked nipples.

A Healthcare Professional or Lactation Professional can give you guidance in breastfeeding your baby and using a breastpump. Their guidance is valuable in terms of the lifetime of health benefits for both you and your baby. A listing of lactation professionals are available at [www.medelalocator.com](http://www.medelalocator.com) or calling 1-800-TELL-YOU (1-800-835-5968).

## 17. Warranty

This product is warranted by Medela, Inc., to the original retail purchaser to be free from defects in material and workmanship for the period of three (3) years from the date of purchase. In the event of a defect or failure to conform to this warranty, Medela will, at Medela's option, repair or replace this product without charge for such replacement, on parts or labor, at the Medela office listed below. The Purchaser shall bear all responsibility and expense for returning this product to Medela, including risk of loss prior to receipt by Medela, and shipping, packing, and insurance costs.

Medela will have no obligation under this warranty to repair any malfunction or damage arising from any misuse, abuse, or alteration of this product. Without limiting the generality of the foregoing, bending, or dropping of this product or its components, visible cracking of the equipment housing will be presumed to be defects resulting from misuse or abuse. Accessories that are not expressly manufactured by Medela and subsequently attached to the equipment will immediately void this warranty.

ANY AND ALL IMPLIED WARRANTIES ARE LIMITED TO THE DURATION OF THREE YEARS FROM DATE OF PURCHASE. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU.

THE LIABILITY OF MEDELA FOR BREACH OF ANY WRITTEN OR IMPLIED WARRANTY IS LIMITED TO REPAIR OR REPLACEMENT OF THIS PRODUCT, AND MEDELA WILL HAVE NO LIABILITY UNDER ANY CIRCUMSTANCES FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATIONS OR EXCLUSIONS MAY NOT APPLY TO YOU.

THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE.

If you wish to make a claim under this warranty, you must return this product to Medela at the address below. You must contact Medela Customer Service and obtain a Return to Merchant Authorization (RMA) Number and reference the number on your return package. The original Invoice Number or Purchase Order Number must be referenced on the RMA. Customer Service toll free number: 1-800-435-8316.

Medela, Inc. – Returns, Door 4501  
1101 Corporate Dr.  
McHenry, IL 60050  
ATTENTION: RETURNS

Call first for authorization number. Returns not accepted without an authorization number.

## 18. Meaning of symbols

The warning symbol identifies all instructions that are important for safety. Failure to observe these instructions can lead to injury or damage to the breastpump. When used in conjunction with the following words, the warning symbol stands for:

### Symbols on the packaging

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates a carton package.

This symbol indicates to keep the device away from sunlight.

This symbol indicates to handle the fragile device with care.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates to keep the device dry.

This symbol indicates that the packaging contains products intended to come in contact with food according to regulation 1935/2004.

This symbol indicates unique GS1 Global Trade Item Numbers (GTIN)

## Symbols on the device

This symbol indicates to follow instructions for use.

This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

This symbol indicates the manufacturer.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This symbol indicates a type BF applied part.

This symbol indicates manufacturer's serial number of the device.

This symbol indicates manufacturer's catalog number of the device.

IP21

This symbol indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water.

This symbol indicates the date of manufacture (four digits for the year and two digits for the month).

This symbol indicates that the Symphony is a class II device.

This symbol indicates fuse (to identify fuse boxes or their location).

This symbol indicates that the fuse box is safety tested.

This symbol indicates open the cover of the mains socket.

Let-down button

On/off button

Indication to press letdown button for Premie+™ program

Vacuum level indicator/Stimulation Phase - in Standard and Premie+ program

Vacuum level indicator/Expression Phase - in Standard program

Vacuum level indication/Expression Phase - in Premie+ program

# 19. EMC/Technical Description

The Symphony® breastpump needs special precautions regarding EMC (Electromagnetic compatibility) and needs to be installed and put into service according to the EMC information provided in these instructions for use.

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations and walkie-talkies can affect the Symphony electric breastpump and should be kept at least a distance 1.0 m away from the equipment.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 1)

Electromagnetic emissions		
The Symphony electric breastpump is intended for use in the electromagnetic environment specified below. The customer or the user of the Symphony electric breastpump should assure that it is used in such an environment.		
Emission tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Symphony electric breastpump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Symphony electric breastpump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

## Warning

**Warning** – The Symphony electric breastpump should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Symphony electric breastpump should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 2)

<b>Electromagnetic immunity</b>			
<p>The Symphony electric breastpump is intended for use in the electromagnetic environment specified below. The customer or the user of the Symphony electric breastpump should assure that it is used in such an environment.</p>			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge ( ESD ) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input/ output lines	Mains Power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line( s ) to line(s) ± 2 kV line( s ) to earth	± 1 kV differential mode ± 2 kV common mode	Mains Power quality should be that of a typical commercial or hospital environment.

## 19. EMC/Technical Description (cont.)

<p>Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p>	<p>&lt;5 % <math>U_T</math> (&gt;95 % dip in <math>U_T</math>) for 0,5 cycle 40 % <math>U_T</math> (60 % dip in <math>U_T</math>) for 5 cycles 70 % <math>U_T</math> (30 % dip in <math>U_T</math>) for 25 cycles &lt;5 % <math>U_T</math> (&gt;95 % dip in <math>U_T</math>) for 5 s</p>	<p>&lt;5 % <math>U_T</math> (&gt;95 % dip in <math>U_T</math>) for 0,5 cycle 40 % <math>U_T</math> (60 % dip in <math>U_T</math>) for 5 cycles 70 % <math>U_T</math> (30 % dip in <math>U_T</math>) for 25 cycles &lt;5 % <math>U_T</math> (&gt;95 % dip in <math>U_T</math>) for 5 s</p>	<p>Mains Power quality should be that of a typical commercial or hospital environment. If the user of the Symphony electric breastpump requires continued operation during power interruptions, it is recommended that the Symphony electric breastpump is powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>3 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>
<p>NOTE <math>U_T</math> is the a.c. voltage prior to application of the test level.</p>			

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 4)

**Electromagnetic immunity**

The Symphony electric breastpump is intended for use in the electromagnetic environment specified below. The customer or the user of the Symphony electric breastpump should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>10 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Symphony electric breastpump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance <math>d = 0,35\sqrt{P}</math></p> <p><math>d = 1,2\sqrt{P}</math> 80 MHz to 800 MHz <math>d = 2,3\sqrt{P}</math> 800 MHz to 2,5 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

## 19. EMC/Technical Description (cont.)

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Symphony electric breastpump is used exceeds the applicable RF compliance level above, the Symphony electric breastpump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Symphony electric breastpump.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 6)

**Recommended separation distances between portable and mobile RF communications equipment and the electric breastpump Symphony**

The Symphony electric breastpump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Symphony electric breastpump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Symphony electric breastpump as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter <i>W</i>	Separation distance according to frequency of transmitter <i>M</i>		
	150 kHz to 80 MHz $d = 0,35\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3\sqrt{P}$
0,01	0,04	0,12	0,23
0,1	0,11	0,37	0,74
1	0,35	1,2	2,3
10	1,11	3,7	7,4
100	3,50	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

# 20. Technical specifications

vacuum (approx.)  
 -50 ... -250mmHg  
 -7 ... -33kPa  
 45 ... 120cpm

100-240V-  
 50/60Hz  
 0.5A


T 0.8A, 250V  
 (slow blow)  
 5 x 20mm L = Low Breaking Capacity

12 V  $\overline{\text{---}}$  Car adapter  
 T 2.5 A  
 (slow blow)

2x6V, 1200mAh  
 Pb (lead acid)  
 Yuasa NP 1.2-6

900 410 410mm

245 215 155mm

3.5kg

2.9kg

Operation

Transport/Storage

Operation

Transport/Storage

Operation/Transport/Storage

# 21. Disposal

## Disposal

---

The unit is made of various metal and plastics. Before disposal, the device is to be rendered unusable and it must not be disposed of as unsorted municipal waste in accordance with local regulations. Use your local return and collection system for waste electrical and electronic equipment (incl. batteries). Improper disposal may have harmful effects on the environment and on public health.

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Printed in the USA. Medela, Symphony and 2-Phase Expression are registered trademarks of Medela Holding AG. Quick Clean, Micro-Steam, Maximum Comfort Vacuum, Preemie and PersonalFit are trademarks of Medela, Inc.
Medela AG
Assembled by Medela, Inc.